# Mattress Care Plus

**06367** Highly elastic quality mattress with great comfort designed for use on high/low beds.

## Core

- The core consists of 1 solid block of 12 cm CMHR foam with a density of 40-45 kg/m³.
- The polyurethane foam is very elastic, fire retardant and thermosetting.
- The mattress retains its elasticity even after frequent use.
- Bacteriostatic treatment against pathogenic and allergenic micro-organisms.
- Smooth and excellent flexibility for use on an adjustable bed.

#### Cover

- The material of the cover is completely adjusted to optimal pressure distribution and a low frictional resistance.
- The cover is waterproof, vapor permeable, breathable and urine resistant.
- The cover is antimicrobial, hypoallergenic and neutralizes pathogenic micro-organisms and allergens.
- Removable, with zipper on a long and short side for easy replacement.
- Washable at 95 °.

-



### **Dimensions**

06367 | Care Plus mattress

- L 195 x W 85 x H 12 cm



# Materials and finishing

- Core: Polyurethane foam CMHR (40-45 kg/m³).
- <u>Cover</u>: Polyurethane coating on a base of polyester
- The mattress and the cover meet fire standard CRIB 5
- The mattress and cover are ecologically produced without CFCs. Neither the production nor the possible destruction processes are environmentally harmful.
- Withstands the common cleaning products.

For more information about the **materials**: please consult our material lists.

