# Talo Classic

12900 The Talo Classic is a comfortable and well-designed relax chair.

The Classic mechanism can be easily adjusted by pushing yourself off a little from the chair, and this without using any levers.

# **Support**

- The support is made of 2 solid beech sides, connected to each other by means of a solid beech cross beam.
- The side parts are connected to each other by means of a glued mortise and tenon joint.
- The sides are firmly connected to each other by means of a solid beech cross beam and a metal frame.
- The ergonomic design of the armrests makes it easy to stand up or sit down.
- Comfortable rounded armrests of  $\pm$  7 cm.
- All edges and corners have been rounded.
- The legs are equipped with synthetic glides.

# Seat

- Removable seat, made of a melaminated MDF board of 12 mm.
- The rounded seat in front stimulates a fluent blood circulation.
- Comfortable seat height of 47 cm.
- Upholstered with fire-retardant foam in 2 different density zones to guarantee an excellent seat comfort.
- The upholstery is finished with imitation leather and decorative stitching.

# **Backrest**

- Made of an ergonomically preformed plywood shell of 12 mm.
- Upholstered with fire-retardant foam in 3 different density zones to guarantee an excellent seat comfort.
- The upholstery is finished with imitation leather and decorative stitching.
- Equipped with an adjustable and removable head/neck cushion.

# Legrest

- Made of a plywood board of 10 mm.
- Upholstered with fire-retardant foam.

# Classic mechanism

- 3 positions: seat position, reading position, and a relax position.
- The relax mechanism first raises the legs, after which it is possible to further incline the entire chair.
- Easy to adjust by pushing yourself off a little.
- Adjustable without needing a lever.
- Guaranteed optimal seat comfort in every position: when the chair is reclined, the angle between the backrest and seat remains constant.
- Optimal leg support, thanks to the telescopic legrest.
- When the chair is put in the seat position, the legrest is turned inwards so that you can stand up without any problem.

### Colours

We would like to refer to our extensive colour gamut. Sample card + samples available on request.



# **Dimensions**

### **12900** | W 68 x D 84 x H 116 cm

- Seat height: 45,5 cm
- Seat depth: 51 cm
- Seat width: 53,5 cmHeight armrests: 64 cm





# Materials and finishing

- Wood: beech (solid), plywood, MDF.
- Metal: epoxy coating (light grey).
- Fire-retardant foam.
- <u>Upholstery</u>: see collection.
- Resistant to the common cleaning products.

For further information about used materials, constructions and maintenance: consult our materialfiles.

# **Options**

- Dinner tray attachment on armrests.
- Seat cushion with visco-elastic foam and a bielastic cover.
- Anatomical reduction cushion.
- Lateral headrest.
- Backrest with removable headrest.
- Integrated auxiliary wheels on the back legs.
- Closed side panels.
- Accessory brackets.
- Restraint bracket.
- Case for IV rod.
- Dinner tray support on backrest.

