# llana Hospital

**12858** The stylish llana bedside table is made of very durable materials such as solid laminate and aluminium.

The bedside table with high drawer can be optionally equipped with a dinner tray.

# Version with high drawer (high drawer, open compartment, drawer)

- The sides are made of 8 mm solid laminate panels that are assembled between 4 aluminium corner profiles.
- Solid laminate drawer fronts.
- Metal design handle.
- An open compartment is accessible from all sides.
- The drawers are made of steel plate and have solid laminated fronts.
- The drawers glide on 2 synthetic rails.
- The upper drawer is equipped with a removable synthetic tray.
- Aluminium and steel parts are protected by an epoxy-coating or an anodisation layer.
- Tabletop with rounded solid laminate corners.
- 4 design castors Ø 65 mm.

# Optional dinner tray Livia

- Folding.
- Easy to incline in both directions.
- The dinner tray has an upright edge.
- The folding part is fixed on a strong aluminium guiding profile that has been finished with an anodisation layer.
- The dinner tray is height-adjustable (between 80 and 110 cm) by means of a gas spring.
- Dimensions of the dinner tray: 63 x 40,5 cm.
- The gas spring is activated by means of a foot pedal.
- The dinner tray is made of fire-retardant ABS.
- The bedside table stands on 5 design castors  $\emptyset$  65 mm.
- The dinner tray can be optionnaly provided with a Click 'n' turn mechanism which allows to move the dinner tray from one side to the other side of the bedside table.

### Handles

We would like to refer to our extensive colour gamut. Sample card + samples available on request.

- Open handle

- Closed handle



### Dimension

12858 | 51 x 42,5 x height 86 cm with dinner tray Livia | 67 x 42,5 x height 86 cm



# Materials and finishing

- Solid laminate.
- Metal: aluminium & steel.
- High-quality self-closing hinges.
- Withstands the common cleaning products.

For more information about the **materials**, **construction and maintenance**: please consult our material lists.

## **Options**

- Livia dinner tray.
- Cylinder lock.
- Electronic lock.
- Gallery.
- Cases for the telephone holder.
- Telephone holder.
- Removable tray with upright edges for the tabletop.
- Bottle tray (for 3 bottles).
- 4 bumper castors.
- 4 castors with brakes.
- Click 'n' turn mechanism.
- USB charger and 2 sockets built into the back.



Colours