



# Handbediening

we create *furniture to heal, to care & to meet*



## Basis functies:

	Back up		Relax position
	Back down		Sleep position
	Knee up		Nurse position <i>only HD</i>
	Knee down		Anti-Trendelenburg
	HL up		UBL
	HL down		Trendelenburg <i>only HD</i>
			Lock function <i>only HD</i>

## Extra functies:

	-L-		Lock All			SAVE Transport position
	-R-		Unlock All (when error = Reset)			Trendelenburg HD handset (magnetic lock)
	-C-		CPR			
	-T-		Transport position			