

# 1. ONE-PIECE SIDEGUARDS

*One-piece sideguards offer maximum safety over three-quarters or the entire bed length.*



## ATMOSPHERE

The **Atmosphere sideguards** offer maximum patient safety in each bed position. The sideguards constitute one unit, protect the patient over 3/4 of the bed length and can be lowered by means of an easy locking mechanism.

The optional **Atmosphere sideguards with 4 siderails** make it possible to use mattresses with a thickness up to 22 cm.



## TRIX

**Three durable aluminium siderails** offer extra safety over the entire bed length. The sideguards can only be lowered by using your both hands. This central safety system prevents the sideguards from being lowered unintentionally.

Max. mattress height Trix sideguards: 16 cm.

### NO RISK OF GETTING STUCK IEC 60601-2-52

All Haelvoet one-piece sideguards comply with safety norm IEC 60601-2:52. This harmonised European norm prescribes the distance between the sideguards and the head section (< 60 mm), between the mattress platform and the sideguards (< 120 mm), between the siderails (< 120 mm), and between the sideguards and the foot section (< 60 mm or > 318 mm).



## TRIX

*telescopic system*

An **integrated telescopic system** makes sure that the Trix sideguards - even when the bed is extended - cover the bed sides over the entire length.

